

SYSTEM OVERVIEW

MFC 1: 10 SLM, OFF ON

Flow: L/min



MFC 2: 10 SLM, OFF ON

Flow: L/min



MFC 3: 30 SCCM, OFF ON

Flow: L/min



MFC 4: 30 SCCM, OFF ON

Flow: L/min



Gas Configuration

Flow Configuration

Calibration Tools

MFC gas

MFC flow

OFF ON



The system overview form is the start up screen for the system. This screen displays the current flow through all mass flow controllers, the fluid flowing through the mfcs, and the type of mfc. The four OFF-ON switches set the flow through each mfc to either zero, or the last input flow (which could be zero). MFC 4 can be a 30 SCCM or 1000 SCCM mfc, which is why the text dictating the type of mfc for MFC 4 is somewhat different than the rest.

The three lower buttons are used to move to gas configuration menu (Form 1), flow configuration menu (Form 2), and calibration tools (Form 4).

GAS CONFIGURATION

MFCs

MFC 1

10 SLM, Dry Air

MFC 2

10 SLM, Humid Air

MFC 3

30 SCCM, Gas

MFC 4

30 SCCM, Gas

Gas

Air(Dry)

Air(Humid)

O2

O3

N2

NO2

CO

CO2

R2

D2

Tank Concentration

PPM

+1000

+100

+10

+1

-1

-10

-100

Reset

Flow Concentration

PPB

+100000

+10000

+1000

+100

+10

-10

-100

-1000

-10000

Reset

Relative Humidity

00.0

%

+10

+1

+0.1

-0.1

-1

-10

Reset

Confirm new settings

Back to main

Concentration

Form 1, the gas configuration menu, setups the inputs into all the mfcs. This menu's information is the basis for most of the system, so it is vital that all pieces are input correctly. A single mfc can be selected at a time, with MFC 1 being locked to dry air and MFC 2 being locked to humid air. MFC 3 and MFC 4 can be changed to whatever gas they are connected to. This setup ensures that users do not try to manipulate the use of MFC 1 or MFC 2, which are used to balance flows after MFC 3 and MFC 4 flows have been considered.

Confirming settings here will not cause any flow to begin, only re balance an already flowing system. Changes in the amount of flow must be done on Form 2, the flow configuration menu.

The bottom right hand corner button sends user back to Form 0. Does not confirm changes.


FLOW CONFIGURATION

Desired Full System Flow

0000 L/min

+1 +0.1 +0.01 -0.01 -0.1 -1 Reset

Individual MFC configuration



System Tempertaure

00 °C

+10 +1 -1 -10 Reset

System Pressure

0000 atm


+1 +0.1 +0.01 -0.01 -0.1 -1 Reset

System Humidity

000 %

+10 +1 +0.1 -0.1 -1 -10 Reset

Back to main

Confirm all changes


Form 2, the system configuration menu, is used to edit the general properties of the system. Temperature and pressure will effect the conversion factors from standard flow to operational flow, as well what temperature the second Arduino will regulate the system to. The full system flow is the combined flow of all four mfcs. System humidity dictates what percentage of the flow should come from MFC 2, which is the humid air flow. Confirming chanes here will result in a change to the system flow, as well as how each mfc calculates how much fluid is flowing through it.

It should be noted that a system flow which causes any mfc to flow less than 1% of its maximum operational flow (without being set to zero), can lead to poorly tuned flows, since the mfcs can not actuate their valves precisely at that level.

The Individual MFC configuration button sends the user to Form 3.

MFC CONFIGURATION

MFCs

| | |
|-------|-------------------|
| MFC 1 | 10 SLM, Dry Air |
| MFC 2 | 10 SLM, Humid Air |
| MFC 3 | 30 SCCM, Gas |
| MFC 4 | 30 SCCM, Gas |

MFC 4 Type

30 SCCM
1000 SCCM

Technical Information

| | | |
|---------------|----------------------|---------|
| Signal In | <input type="text"/> | 0-65535 |
| Signal Out | <input type="text"/> | 0-65535 |
| Standard Flow | <input type="text"/> | L/min |

Collect Tech Info

MFC Flow

0000

L/min

Signal Offset

00000

0-65535

+100 +10 +1 -1 -10 -100 Reset

Signal Gain

0000

per Volt

+100 +10 +1 -1 -10 -100 Reset

Confirm new settings

Back to main

Form 3, the mfc configuration menu, allows users to edit individual mfc settings. Technical Information for an mfc can be found here as well, though is unreliable as of now. Standard flow is the individual mfc's flow if it was at 1 atm and 25 C with the same signal strength. MFC 4's type can be changed from this menu too (30 SCCM to 1000 SCCM or vice versa). Confirming changes on this menu will change the signal strength sent to the affected mfc.

| Concentration MFC 3 | | Concentration MFC 4 | | Increment Time | |
|-----------------------|--------------------------|-----------------------|--------------------------|------------------------------|----------------------|
| OFF | | OFF | | 5 minutes at each increment | |
| Low | <input type="text"/> PPB | Low | <input type="text"/> PPB | 15 minutes at each increment | <input type="text"/> |
| Medium | <input type="text"/> PPB | Medium | <input type="text"/> PPB | 40 minutes at each increment | |
| High | <input type="text"/> PPB | High | <input type="text"/> PPB | | |
| Increment through all | | Increment through all | | Purge Between Increments | |
| Custom Incrementation | | Custom Incrementation | | OFF | |
| | | | | ON, 1 minute | |
| | | | | ON, 5 minute | |
| | | | | Time Elapsed | |
| | | | | <input type="text"/> | |
| | | | | Confirm and Begin | |
| | | | | Stop | |
| | | | | PAMM ID's | |

| Relative Humidity | | Temperature | | In Progress Stats | |
|-----------------------|-----|-----------------------|--|-------------------|--------------------------|
| 15 % | Off | 5 C | | MFC 3 | <input type="text"/> PPB |
| 40 % | | 20 C | | MFC 4 | <input type="text"/> PPB |
| 65 % | | 35 C | | RHum | <input type="text"/> % |
| Increment through all | | Increment through all | | Temp | <input type="text"/> C |

| Concentration of fluid at each specified level. | |
|---|--|
| <input type="text"/> | Information for the ongoing step of the calibration. |
| <input type="text"/> | Time since calibration began. |

Form 4, the calibration setup menu, is rather robust. A set of predetermined concentrations for MFC 3 and MFC 4 can be selected, or users can go to Form 5, the custom increment menu, to setup their own. Relative humidity as well as temperature refer to the system's properties, not any individual mfc. The length of each step of the calibration is selected under Increment Time.

Purge sets all flows besides dry air to zero. Thus "flushing" all mixtures besides air out of the system. This is good to use in order to see a stark contrast between each step, but causes steps to take longer to reach their stable value.

Confirm and Begin starts the calibration. To ensure no calibration stacking, the only responsive button when a call is in progress is the Stop button. PAMM ID's sends you to Form 6, PAMM ID input.

MFC 3

Initial Concentration

PPB

| | | | | | | | | | |
|---------|--------|-------|------|-----|-----|------|-------|--------|-------|
| +100000 | +10000 | +1000 | +100 | +10 | -10 | -100 | -1000 | -10000 | Reset |
|---------|--------|-------|------|-----|-----|------|-------|--------|-------|

Final Concentration

PPB

| | | | | | | | | | |
|---------|--------|-------|------|-----|-----|------|-------|--------|-------|
| +100000 | +10000 | +1000 | +100 | +10 | -10 | -100 | -1000 | -10000 | Reset |
|---------|--------|-------|------|-----|-----|------|-------|--------|-------|

Increment

PPB

| | | | | | | | | | |
|---------|--------|-------|------|-----|-----|------|-------|--------|-------|
| +100000 | +10000 | +1000 | +100 | +10 | -10 | -100 | -1000 | -10000 | Reset |
|---------|--------|-------|------|-----|-----|------|-------|--------|-------|

Mirror
Incrementation

Confirm Custom
Increments

Return to
Calibration Setup

MFC custom
increment is being
made for.

Concentration for
desired property.

Form 5, the custom incrementation menu, allows users to select what concentrations the calibration will step through for a given mfc. The Increment field should be perfectly divisible by Final Concentration minus Initial Concentration, otherwise the final concentration will be over shot on the last step; therefore it will be skipped. Mirror Incrementation can be selected to have the calibration step up through these custom increments, then back down; instead of ending when the final concentration is reached and its increment time finished. Return to Calibration Setup does not save any changes made, sends user back to Form 4. Confirm Custom Increments saves changes, and sends user back to Form 4.

The final concentration must be higher than the initial concentration before the initial can be increased.

Number of PAMM IDs

00

PAMM IDs in comma seperated string

| | | |
|--|--|--|
| | | |
| | | |
| | | |
| | | |
| | | |

Delete last ID

| | | |
|---------|---|----|
| 7 | 8 | 9 |
| 4 | 5 | 6 |
| 1 | 2 | 3 |
| Confirm | 0 | <- |

Back to Calibration Setup

List of all PAMM ID's input.

Form 6, the PAMM ID input screen, allows users to input which PAMMs will be used for the upcoming calibration. The PAMM ID's will show up as a comma separated list. The total number of input ID's is shown in the upper left hand side. Deleting a PAMM ID is done by pressing Delete last ID button. A PAMM ID can not be deleted from the middle of the list.

This is the only menu where the bottom right hand side circular button does not send you back to Form 0. Instead it sends you back to Form 4. All changes must be confirmed before pressing this button, if the user desires them to take effect.